



Adam Equipment's MDW mechanical scales offer an affordable solution for checking the weight of fitness center or gym members. Wheels allow easy mobility and a notched counter weight clearly shows the numbers. MDW scales feature rugged construction with a metal base and a non-slip weighing platform to ensure stability and comfort. Models include an adjustable height rod to simultaneously measure height and weight.

Affordable, easy-to-operate scale measures height and weight



Features

- Height measurements from 60cm to 212cm
- Integrated wheels permit easy movement
- Anti-slip platform ensures safety when stepping on and off the scale
- Height rod facilitates measuring the height of children and adults while weighing
- Metal construction for base and pillar
- Zero adjustment

Applications

- Weighing
- Height measurement

About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Model	MDW 160M	MDW 200M
Capacity	160kg	200kg
Readability	0.1kg	
Repeatability (S.D.)	0.2kg	
Linearity (+/-)	0.2kg	
Pan Size	375x275mm	
Weighing Units	kg, lb	
Stabilisation Time (sec)	3	
Body Mass Index [BMI]	No	
Height Rod	Yes	Yes
Operating Temperature	5° to 40°C	
Overall Dim.	530x275x1485mm (wxdxh)	
Net Weight	15kg	

Supplied by